

BLS ACCESSORIES

C90 – WALL MOUNTED CONTAINER FOR FULL FACE MASKS

The C90 wall container allows to store a full-face mask BLS 5000 / 3000 series and its filters in a quick and functional way. It is equipped with three rear holes for wall mounting, vertical opening and interlocking closure. Ideal in situations where it is necessary to be able to readily dispose of the protective device in a specific place.



MATERIALS

CONTAINER	Polystyrene (PS)
COVER	Polyvinyl chloride (PVC)
HINGES	Steel

TECHNICAL DATA

CODE	PRODUCT	COMPATIBLE MASK	WEIGHT (g)
8001027	C90	BLS 5000 / 3000 SERIES	75

IMPORTANT

BLS declines any responsibility, direct or indirect, from any misuse of both devices and instructions. User is responsible for the determination of product compliance with the intended use.

BLS ACCESSORIES

K15 / K16 / K19 – PROTECTIVE FILMS FOR FULL FACE MASKS



K15

Protective film in polycarbonate for the visor of BLS 5000 series full face masks. Its side holes snaps on the side hooks of the mask. Protects from mechanical wear, in case of sandblasting or generally in presence of abrasive substances.

K19 and K16

Protective film in cellulose acetate for the visor of BLS 5000 series full face masks (K19) or BLS 3000 / 2000 series full face masks (K16). It is applied thanks to adhesive side strips and it has a flap for easy removal. Protects from chemical wear, in contexts such as industrial painting.



TECHNICAL DATA

CODE	PRODUCT	COMPATIBLE MASK	MATERIAL	THICKNESS (µm)	QUANTITY/PACK
8001026	K15	BLS 5000 SERIES	Polycarbonate (PC)	750	1
8001112	K16	BLS 3000 / 2000 SERIES	Cellulose acetate	120	5
8001125	K19	BLS 5000 SERIES	Cellulose acetate	120	5

IMPORTANT

BLS declines any responsibility, direct or indirect, from any misuse of both devices and instructions. User is responsible for the determination of product compliance with the intended use.

BLS ACCESSORIES

C10 – TEXTILE HARNESS FOR FULL FACE MASKS

This textile harness can replace the standard rubber one making wearing operations more comfortable, especially in case of long hair.



MATERIALS

LACES Polypropylene (PP)
NET Polypropylene (PP)

TECHNICAL DATA

CODE	PRODUCT	COMPATIBLE MASK	WEIGHT (g)
8001006	C10	BLS 5000 / 3000 SERIES	60

IMPORTANTE

BLS declines any responsibility, direct or indirect, from any misuse of both devices and instructions. User is responsible for the determination of product compliance with the intended use.

BLS ACCESSORIES

C22 – FRAME FOR PRESCRIPTION LENSES FOR FULL FACE MASKS

Frame for prescription lenses for BLS 5000 / 2000 series full face masks.



MATERIALS

FRONT Polycarbonate (PC)
SIDES ABS + Polycarbonate (PC)

TECHNICAL DATA

CODE	PRODUCT	COMPATIBLE MASK	WEIGHT (g)	VISUAL FIELD	COLOUR
8001060	C22	BLS 5000 / 2000 SERIES	17	> 40	BLACK

IMPORTANT

BLS declines any responsibility, direct or indirect, from any misuse of both devices and instructions. User is responsible for the determination of product compliance with the intended use.

BLS ACCESSORIES

C41/C43 – BAGS FOR HALF MASKS AND FULL FACE MASKS



C41

This bag allows the full face mask and its filters to be stored at the end of the shift. Equipped with zip closure and shoulder strap, it is easily washable.

C43

This bag allows the half mask and its filters to be stored at the end of the shift. Equipped with zip closure, it is easily washable.



MATERIALS C41

COATING	Polyvinyl chloride (PVC)
PADDING	Polyethylene (PE)
SHOULDER STRAP	Polypropylene (PP)

MATERIALS C43

COATING	Nylon
SHOULDER STRAP	Polypropylene (PP)

TECHNICAL DATA

CODE	PRODUCT	TYPE OF MASK	COMPATIBLE MASK	COLOUR	WEIGHT (g)
8001017	C41	FULL FACE MASKS	BLS 5000 / 3000 / 2000 SERIES	BLACK + BLS LOGO	194
8002122	C43	HALF MASKS	BLS 4000next / SGE46 SERIES	GREY + BLS LOGO	65

IMPORTANT

BLS declines any responsibility, direct or indirect, from any misuse of both devices and instructions. User is responsible for the determination of product compliance with the intended use.

BLS ACCESSORIES

Q530/Q409 – DISINFECTANTS FOR HALF MASKS AND FULL FACE MASKS



Q530 - DISINFECTANT WIPES

Wipes, in scented non-woven fabric, suitable for the rapid and effective cleaning and disinfection of plastic surfaces of half masks and full face masks.

COMPOSITION: Benzalkonium chloride, O-phenylphenol, Alkylisoquinoline bromide, chamomile extract in alcoholic solution, detergent excipient, purified water.

Q409 - 150 ml DISINFECTANT SPRAY BOTTLE

Quickly disinfects and deodorizes devices. Active on the Coronavirus family. Tests performed by an accredited EU laboratory. Do not use on facepieces. It does not perform a cleaning action.

COMPOSITION: Benzalkonium chloride, O-phenylphenol, isopropyl alcohol, perfume, propellants.



TECHNICAL DATA

CODE	PRODUCT	MOQ	MEDICAL AND SURGICAL UNIT
8013350	Q530	Box of 12 packs (10 tissues each)	Ministry of Health Registration 12-671
8013351	Q409	1 spray bottle	Ministry of Health Registration 13-431

IMPORTANT

BLS declines any responsibility, direct or indirect, from any misuse of both devices and instructions. User is responsible for the determination of product compliance with the intended use.